

MEDICAM

Video Laryngoscope SC-01



MEDICAM'S Feature Rich Video Laryngoscope SC01

Intubation procedure require skills and secondary confirmation device to ensure correct positioning of the tube. MEDICAM Video laryngoscope is a dependable device that helps anesthesiologist to carry out the complex procedure with ease.

The product looks like a normal laryngoscope retaining all the advantages. But what makes it class apart is the small camera fitted on it. One tiny addition that makes a big difference.

Features

- Waterproof camera
- CMOS Sensor
- Inbuilt LED Illumination
- Video recording/Image storing
- Instant playback
- Easy data storage on MicroSD card
- 3.5 inch multifunction screen
- Specially designed concave blade
- 180 ° rotating LCD screen
- Multiple blades option of different types & sizes
- Downloading option
- AV out on TV/AV monitor
- Strong perfect size handle with coating
- Steel blades
- One year warranty
- Neonatal, Pediatric and Adult Blades sets available



Camera

Probe Dimension	Ø5.5mm
Imaging Sensor	1/12" VGA sensor
Effective Pixels	720(H) x 525(V)
Focus Distance	30mm-80mm
View Angle	60 degree F.O.V.
Working Voltage	5V
Working Current	600MA/H

Voltage Adapter

Input Voltage	110-220V
Output Voltage	5V/2A

Screen Specifications

Screen Battery	Built-in Li-battery 1200MAH 3.7V
Screen Working Time	1-2 Hour
Screen Charging time	4-6 Hours
Screen	3.5" TFT LCD
Image Format	JPEG
Video Format	AVI
Functions	Real time clock with photo and video. Time label Image saving Video Recording Magnification Automatic recognition of NTSC/PAL Composite video output
Memory Card	Capacity of max. 32G mini SD/TF card
USB Interface	Downloading stored data
Data Output	Video output, stored data downloading on computer
Language	English
Host Dimension	Screen : 112 X 82 X 25.5mm
Battery	Li Ion



An ISO 9001:2015 Company
An ISO 13485:2012 Company

Shop No. 20/22, Bhivandiwala Bldg, Dr. Ambedkar Road, Kalachowki, Mumbai - 400033
Email: info@hemantmedicam.com, sales@medicam.co.in, hemantmedicam@yahoo.co.in
Ph No.: +91 22 24700872, 8767303439
Website: www.medicam.co.in