


USB Video Proctoscope

1 Image = 1000 Words



Invest in Tomorrow, Practice it Today.

Features

- Small and Lightweight
- No fiber optic and light source
- Portable
- Works almost with all proctoscopes
- Built in adjustable Illumination
- Convenient to use
- Integrated software to freeze and store
- USB connectivity
- No power supply or batteries required
- Auto Focus
- Anti-Fog
- Waterproof
- Shock Proof
- Image Saving
- Slider to adjust the light illumination
- Capturing buttons (for capturing image)
- USB cable (1.5 meters)

Specification

- Sensor type - CMOS sensor
- Lens - Pin hole lens -(Antifog)
- Illumination - 6 White LEDs
- Power Source - USB
- Focus Range - Auto Focus (1.5-3 inches)
- Integrated Software
- Software compatibility = Win XP, Win7, Win 8 & Win 10
- Storage = on Hard disk
- IP67 level of Water Proof.
- Camera probe diameter= 5.5 mm Φ
- Imaging Sensor = 1/12" VGA sensor
- Effective Pixels = 720(H) X 525(V)
- View Angle = 60° FOV
- Working Voltage = 5V
- Working Current = 600 mA
- Image format = JPEG(Real Time Image.)
- 6 Led (With Intensity Adjustment)
- Length of Camera Probe (Metal Flexible probe=8.5 inches)

ADULT COUPLERS DIMENSIONS

Height	45.0 mm
Lower width	32 mm
Upper width	23 mm
Diameter	6.5 mm

PAEDIATRIC COUPLERS DIMENSIONS

Height	41 mm
Lower width	30.5 mm
Upper width	18 mm
Diameter	6.5 mm

MEDICAM

SH. NO. 20/22, BHIWANDIWALABLDG, DR. AMBEDKARMARG, CHINCHPOKALI (EAST), KALACHOWKI, MUMBAI 400033.

Website: www.medicam.in • Customer Care No.: 022-24700872

• Email: support@medicam.in/info@medicam.in/info@hemantmedicam.com