


# VIDEO STYLET

**An Innovative Advancement in Intubation** The Video Stylet by MEDICAM is a groundbreaking device designed to significantly enhance the intubation process. Unlike a conventional stylet, this innovative tool offers multiple advantages that streamline clinical procedures and improve patient outcomes.

## Key Features :

- **Integrated Micro Camera:** A compact, waterproof, and anti-fog micro camera is seamlessly embedded into the stylet, providing real-time visuals that aid in precise intubation, especially in challenging airway scenarios.
- **Durable Stainless Steel Body:** Constructed with high-quality stainless steel, the Video Stylet is not only strong and durable but also easy to clean and sterilize—ensuring safety and longevity.
- **Anatomically Optimized Curve:** The stylet's ergonomic curve is specifically designed to match the natural oral anatomy, allowing for smoother, quicker, and less traumatic intubation.
- **Ultra-Slim Design:** Its remarkably small diameter makes it suitable for use in most difficult airway cases, expanding its versatility across a wide range of patients and clinical situations.

The Video Stylet is the perfect combination of precision, durability, and innovation—setting a new standard in airway management.



## **F** EATURES

- Anti fog
- Perfect shaped curve to suit oral cavity
- Waterproof camera
- Strong steel body
- User-friendly
- In-built Illumination
- Useful in difficult intubations
- Video recording



## **S** pecifications

- Video Stylet Specification.
- Probe Dimension- 5 mm OD
- Imaging sensor-CMOS
- View angle- 88°
- Working Voltage - 5v
- Connectivity type -USB
- Weight of product -280gm.
- Outer Material of product- SS304
- Water resistance- IP 67
- ProductDimension-length140cm
- Output Resolution= 1280\*720 pixels

